

# 4Healthcare

## Instrument Dressing Trolleys



### Assembly & Care Instructions

1. Remove handle screws and replace correct way round as shown.  
Do not over tighten.
2. Check castors roll freely and are lint free.
3. Clean with a mild detergent or a diluted solutions of sodium hypochlorite using a soft cloth or sponge.
4. Maximum load capacity:  
Trolley - 20kgs (10kg/shelf)  
Drawer - 4kgs



4/10 Uppill Place Wangara WA 6065  
Freecall: 1300 525 520 Fax: (08) 9309 2953  
Website: [www.4healthcare.com.au](http://www.4healthcare.com.au)  
Email: [sales@4healthcare.com.au](mailto:sales@4healthcare.com.au)

\*Due to 4Healthcare's commitment to continuous improvement, our products are subject to change without notice.  
Form #: 4HC0085